product preview ws3120



CHALLENGES

- Protective aprons are damaged as a result of improper storage.
- Wall mounted storage systems that pull out of the wall after minimal use.
- Inadequate mounting footprint and mounting hardware.
- Hard to clean metal spot welds, chipping paint that exposes bare metal where bacteria can grow.
- Unattractive design limits mounting locations when a warm / comfortable / visual environment is desired.
- · Flimsy hangers that are awkward to access and return the protective garment.

SOLUTION

• wall**shoulders**[™] by DeBroeck Company; with our unique patent pending design and DuPont[™] Corian[®] raw material, we are proud to bring you the next generation of x-ray garment storage solutions for the healthcare industry.

BENEFITS

- Our ws3120 provides an attractive and durable storage for your back Velcro poncho style panoramic x-ray apron.
- Our unique design provides easy access to remove or replace the apron from its protective and elegant storage
- The full top radius provides full contact support protecting the apron from becoming damaged as a result of improper storage.
- By adding our gs1126 your technician will have a convenient storage for their x-ray protective gloves. (see gs1126)
- The overall width of 20 1/2 inches makes finding an installation location a snap.
- wall**shoulders**[™] are inviting to the touch as well as elegant: The four colors available will accent any rooms design.
- All garment contact points are set away from the mounting surface, this feature keeps the garment from rubbing the mounting surface while in storage.
- Every wallshoulders[™] is manufactured using the highest quality **DuPont[™] Corian**® raw material, the same material used in healthcare facilities around the world. **DuPont[™] Corian**® is non pours and is resistant to bacteria and staining, it is ideal for use in a bacteria intense environment.



DuPontTM Corian® is...

the perfect surface for healthcare

• Non-Porous

Will not support the growth of mold, mildew or bacteria.

• Easy To Clean

No spot welds, creases or gaps; maintenance cost and time are greatly reduced.

• Stain Resistant

Stains cant permeate and are easily removed.

• Durable

Tough and long lasting, stands up to heavy use.

Renewable

Marks and scratches can be removed easily to restore it to its original condition.

Warm And Inviting to the touch

The practical surface material that creates a unique, beautiful statement.

Code Compliant

Class I(A) fire rated and NSF/ANSI 51 Certified for food contact of all food types. Can be molded and formed with no sharp edges or corners.

How Do You Protect Your X-Ray Aprons?

A secure Installation is a priority to us

By following the clear and precise Installation Instructions provided, Your new wall**shoulders**™ will stay put!

No More Wobbling / Hanger / Headaches.

Please copy and paste the following link to your browser for a short animated installation demonstration of the TOGGLER® hollow wall anchor provided in every wallshoulders[™] installation packet.

http://www.toggler.com/products_hwm_installation.html

While there click on the technical information tab for hold ratings



Other Corian®-Based Healthcare Products By DeBroeck™ Company

- gs1126 X-Ray Protective Glove Storage, accepts up to a 15 inch pair of gloves.
- ws3130 Panoramic Shoulder Velcro X-Ray Apron Storage.
- ws2120 Loop or Grommet Style X-Ray Apron Storage.
- Special Order Custom Solid Surface Counter tops, call or e-mail us for details.
- Look for updates on our new products for 2008.

SAMPLES / LITERATURE / QUESTIONS / SPECIAL REQUESTS

To request Actual Color / Material samples and if you have any Questions or Special Requests, please visit us at: www.wallshoulders.com

Please visit our FAQ's section, while there you can submit your Question or Browse the topics already submitted. If you prefer you can also e-mail us direct at: info@wallshoulders.com * or call us at: (509) 525-1349

