

product preview ws2120



wall**Shoulders**[®]
the solution to your x-ray garment storage

CHALLENGES

- Protective aprons are damaged as a result of improper storage.
- Wall mounted storage systems that pull out of the wall after minimal use.
- Inadequate mounting footprint and mounting hardware.
- Hard to clean metal spot welds, chipping paint that exposes bare metal where bacteria can grow.
- Unattractive design limits mounting locations when a warm / comfortable / visual environment is desired.
- Flimsy hangers that are awkward to access and return the protective garment.

SOLUTION

- wall**Shoulders**[®] by DeBroeck Company; with our unique patent pending design and DuPont[™] Corian[®] raw material, we are proud to bring you the next generation of x-ray garment storage solutions for the healthcare industry.

BENEFITS

- Our ws2120 provides an attractive and durable storage for your grommet or loop style x-ray apron.
- Our unique design provides easy access to remove or replace the apron from its protective and elegant storage with ease.
- By adding our gs1126 your technician will have a convenient storage for their x-ray protective gloves. (see gs1126)
- Shoulder wings support the aprons shoulders from being scrunched up or creased while not in use.
- The overall width of 20 1/2 inches makes finding an installation location a snap.
- wall**Shoulders**[®] are inviting to the touch as well as elegant: *The four colors available will accent any rooms design.*
- All garment contact points are set away from the mounting surface, this feature keeps the garment from rubbing the mounting surface while in storage.
- Every wall**Shoulders**[®] is manufactured using the highest quality **DuPont[™] Corian[®]** raw material, the same material used in healthcare facilities around the world. **DuPont[™] Corian[®]** is non porous and is resistant to bacteria and staining, it is ideal for use in a bacteria intense environment.



DuPont™ Corian® is...
the perfect surface for healthcare



2
0
0
8

• **Non-Porous**

Will not support the growth of mold, mildew or bacteria.

• **Easy To Clean**

No spot welds, creases or gaps; maintenance cost and time are greatly reduced.

• **Stain Resistant**

Stains cant permeate and are easily removed.

• **Durable**

Tough and long lasting, stands up to heavy use.

• **Renewable**

Marks and scratches can be removed easily to restore it to its original condition.

• **Warm And Inviting to the touch**

The practical surface material that creates a unique, beautiful statement.

• **Code Compliant**

Class I(A) fire rated and NSF/ANSI 51 Certified for food contact of all food types. Can be molded and formed with no sharp edges or corners.

A secure Installation is a priority to us

By following the clear and precise Installation Instructions provided,
Your new wall**Shoulders**® will stay put!


No More Wobbling / Hanger / Headaches.

Please copy and paste the following link to your browser for a short animated installation demonstration of the TOGGLER® hollow wall anchor provided in every wall**Shoulders**® installation packet.

http://www.togglor.com/products_hwm_installation.html

While there click on the technical information tab for hold ratings

How Do You Protect Your X-Ray Aprons?



Other Corian®-Based Healthcare Products By DeBroeck™ Company

- gs1126 X-Ray Protective Glove Storage, accepts up to a 15 inch pair of gloves.
- ws3120 Panoramic Poncho X-Ray Apron Storage.
- ws3130 Panoramic Shoulder Velcro Style X-Ray Apron Storage.
- Special Order Custom Solid Surface Counter tops, call or e-mail us for details.
- Look for updates on our new products for 2008.

SAMPLES / LITERATURE / QUESTIONS / SPECIAL REQUESTS

To request *Actual Color / Material samples* and if you have any **Questions or Special Requests**, please visit us at:
www.wallshoulders.com

Please visit our **FAQ's** section, while there you can submit your **Question** or **Browse** the topics already submitted.
If you prefer you can also e-mail us direct at: info@wallshoulders.com * or call us at: (509) 525-1349