

product preview gs I 126



wall**Shoulders**®
the solution to your x-ray garment storage

CHALLENGES

- A storage system that can store one pair of x-ray protective gloves (gs I 126) and one protective x-ray garment (ws3 I 30) while providing full shoulder support.
- Reducing the stress and breakdown of the slanted style apron shoulders when stored on straight bar style racks.
- Wall mounted storage systems that pull out of the wall after minimal use.
- Inadequate mounting footprint and mounting hardware.
- Hard to clean metal spot welds, chipping paint that exposes bare metal where bacteria can grow.
- Un-attractive design limits mounting locations when a warm / comfortable / visual environment is desired.
- Flimsy hangers that are awkward to access and return the protective garment.
- A glove storage (gs I 126) that can be installed piggyback with the protective apron storage (ws3 I 30) or installed separately as a glove only storage.

SOLUTION

- wall**Shoulders**® by DeBroeck Company; with our unique patent pending design and DuPont™ Corian® raw material, we are proud to bring you the next generation of x-ray garment storage solutions for the healthcare industry.

BENEFITS

- Our gs I 126 provides easy access to retrieve or return protective gloves even while the protective apron is in storage.
- wall**Shoulders**® gs I 126 is designed specifically for the storage of *x-ray protective gloves*. It is recommended to install the gs I 126 piggyback to our ws3 I 30, providing a complete storage system for 1 X-ray Apron & 1 Pair of Protective Gloves.
- The gs I 126 can also be installed separately in a different location if required or preferred
- The four DuPont™ Corian® colors available will accent any rooms design.
- Every wall**Shoulders**® is manufactured using the highest quality **DuPont™ Corian®** raw material, the same material used in healthcare facilities around the world.
- **DuPont™ Corian®** is non porous and is resistant to bacteria and staining, it is ideal for use in a bacteria intense environment.
- wall**Shoulders**® gs I 126 can also be installed with any of our other wall**Shoulders**® series garment hangers.



DuPont™ Corian® is...
the perfect surface for healthcare



2
0
0
8

• **Non-Porous**

Will not support the growth of mold, mildew or bacteria.

• **Easy To Clean**

No spot welds, creases or gaps; maintenance cost and time are greatly reduced.

• **Stain Resistant**

Stains can't permeate and are easily removed.

• **Durable**

Tough and long lasting, stands up to heavy use.

• **Renewable**

Marks and scratches can be removed easily to restore it to its original condition.

• **Warm And Inviting to the touch**

The practical surface material that creates a unique, beautiful statement.

• **Code Compliant**

Class I(A) fire rated and NSF/ANSI 51 Certified for food contact of all food types. Can be molded and formed with no sharp edges or corners.

A secure Installation is a priority to us

By following the clear and precise Installation Instructions provided, Your new wall**Shoulders**® will stay put!

No More Wobbling / Hanger / Headaches.

Please copy and paste the following link to your browser for a short animated installation demonstration of the TOGGLER® hollow wall anchor provided in every wall**Shoulders**® installation packet.

http://www.toggler.com/products_hwm_installation.html

While there click on the technical information tab for hold ratings

How Do You Protect Your X-Ray Aprons?



Other Corian®-Based Healthcare Products By DeBroeck™ Company

- ws3130 Panoramic Shoulder Velcro X-Ray Apron Storage.
- ws3120 Panoramic Poncho X-Ray Apron Storage.
- ws2120 Loop or Grommet Style X-Ray Apron Storage.
- Special Order Custom Solid Surface Counter tops, call or e-mail us for details.
- Look for updates on our new products for 2008.

SAMPLES / LITERATURE / QUESTIONS / SPECIAL REQUESTS

To request *Actual Color / Material samples* and if you have any **Questions or Special Requests**, please visit us at:
www.wallshoulders.com

Please visit our **FAQ's** section, while there you can submit your **Question** or **Browse** the topics already submitted.
If you prefer you can also e-mail us direct at: info@wallshoulders.com * or call us at: (509) 525-1349